

MDPLANNER®

PRODUCTS OVERVIEW



Nutricare

PRODUCTS LINE-UP

BEAUTY DEVICE



Nutricare

MDPLANNER I3

BEAUTY DEVICE



KEY FUNCTION

Galvanic Ion Treatment

It helps active ingredients from skincare cosmetics penetrating into skin layer.

LED Light Therapy (3 Colors)

- (1) Blue LED on Cleansing Mode
- (2) Green LED on Nutrition Mode
- (3) Red LED on Lifting Mode



Vibration Massage

It promote the absorption of active ingredients and help for the better blood circulation which makes skin tone healthier.

COMPOSITION

Main Unit, AAA Battery*2ea

**Nutricare**

MDPLANNER I7

BEAUTY DEVICE



KEY FUNCTION

Galvanic Ion Treatment

It helps active ingredients from skincare cosmetics penetrating into skin layer.

LED Light Therapy (4 Colors)

- (1) Blue LED on Cleansing Mode
- (2) Green LED on Nutrition Mode
- (3) Red LED on Lifting Mode
- (4) Yellow LED on Brightening Mode



Heat-Up Thermotherapy

On thermotherapy mode, it activates cells and promote blood circulation.

Vibration Massage

It promotes the absorption of active ingredients and help for the better blood circulation which makes skin tone healthier.



COMPOSITION

- Main Unit, Charging Cradle, Adaptor, USB Cable, Package Box, Ion Band(Optional for Professional Use)

MDPLANNER E3

BEAUTY DEVICE



KEY FUNCTION

4-Minute Anti-Wrinkle Program

EMS (3 minutes) + Galvanic Ion Therapy (1 minute) = 4 minutes

EMS (Electrical Muscle Stimulation)


- We can specially care muscles on eye zone such as we work out and exercise to firm up our body muscle.
- Customized low-frequencies stimulate the muscle layer to tighten and lift up like 20s skin condition.
- Wrinkles and Drooping eye-bags can be removed
- 4 kinds of Irregular low-frequency stimulates the muscle without muscle adaptation problem.

Galvanic Ion Treatment

It helps active ingredients from skincare cosmetics penetrating into skin layer.

COMPOSITION

Main Unit, AA Battery*1ea, User Manual, Package Box



MDPLANNER E5

BEAUTY DEVICE



KEY FUNCTION

4-Minute Anti-Wrinkle Program

- Intensive EMS Treatment for 4 minutes + Red LED Therapy
- Ergonomic Body and 2 heads are perfectly massage whole face.

EMS (Electrical Muscle Stimulation)


- We can manage the face muscle such as we work out and exercise to firm up our body muscle.
- Customized low-frequencies stimulate the muscle layer to tighten and lift up like 20s skin condition.
- Wrinkles and Drooping eye-bags can be removed
- Irregular low-frequency stimulates the muscle without muscle adaptation problem.

< Wave of Irregular Low-Frequency >



COMPOSITION

- Main Unit, USB Charging Cable, User Manual, Package Box




MDPLANNER R3

BEAUTY DEVICE



KEY FUNCTION

- **Total Anti-aging Care System (4 minutes)**
RF(High Frequency) + Galvanic Ion Treatment + Red LED Therapy
- **Bi-polar RF – high technology of skin treatment.**
By applying Bi-polar High Frequency to generate the heat on the deep skin layer (not skin surface), it promotes the synthesis of skin protein, stimulate the blood circulation, and help to grow collagen and elastin to realize the younger skin age ultimately.
- **Galvanic Ion Treatment**
It helps active ingredients from skincare cosmetics penetrating into skin layer, making your skin healthier and brighter.



COMPOSITION

- Main Unit, USB Charging Cable, User Manual, Package Box



MDPLANNER R7



BEAUTY DEVICE



KEY FUNCTION

The Most Powerful Facial Anti-aging Device

- World's Best Innovative Skincare Technology
- Intensive RF(High Frequency) Treatment + Galvanic Ion Treatment + Near-infrared Ray Therapy + 3 Colors of LED Therapy (Blue, Yellow, Red)

Bi-polar RF(High Frequency)


By applying Bi-polar High Frequency to generate the heat on the deep skin layer (not skin surface), it promotes the synthesis of skin protein, stimulate the blood circulation, and help to grow collagen and elastin to realize the younger skin age ultimately.

Near-Infrared Ray Therapy

NIR is penetrated into deep skin to activate the skin cell and blood circulation to make the skin healthy.

Galvanic Ion Treatment

It helps active ingredients from skincare cosmetics penetrating into skin layer, making your skin healthier and brighter.



COMPOSITION

- Main Unit, Charging Cradle, Adaptor, USB Cable, Package Box, Ion Band(Optional for Professional Use)


< High Resolution LED Colors >



Red LED

Yellow LED

Blue LED



MDPLANNER C4

BEAUTY DEVICE




KEY FUNCTION

Patented 4D Electric Cleaner (Patent number:10-1663212)

4D cleansing brush performs in 4 directions (right-and-left and up-and-down) for the most effective deep cleansing in pores.


Finest Filaments in the World!

MDPLANNER C4 adopts the world's finest filaments to reduce any friction or skin irritation.



Easy to Wash Out

The brush head is detachable for easier washing-out after facial cleansing.



COMPOSITION

Main Unit, AAA Battery*2ea

**Nutricare**

MDPLANNER B9

BEAUTY DEVICE



KEY FUNCTION

- IV+ EMS + eV + Red LED Light Therapy**
Synergy effect by perfect Combination of 4-kind body fitness
- IV (Interferential Vibration)**
It relaxes muscles before intensive EMS treatment and boost the slimming effect by warming up.
- EMS (Electrical Muscle Stimulation)**
5 different frequencies strongly build up the muscle layer and burn up fat and calories.
- eV (Electrical Vibration)**
After EMS treatment, eV electrical massage system takes care of stimulated muscles by using 4 different professional hand massage technique.

COMPOSITION

- Main Unit, Cradle, USB Cable, User Manual, Package Box**



MDPLANNER H2-1000

BEAUTY DEVICE



KEY FUNCTION

By drinking H₂ Water, Gain your Youth and Beauty back!

- Hydrogen Water neutralizes Hydroxyl free radicals, the most dangerous free radicals, which directly attack and damage cell, protein, DNA and others.



Extraordinary Hydrogen Water Maker!

- 1,000ppb of hydrogen intensity within 5 minutes
- [Patent #. 10-1844968] New type of twist-and-open lid
- BPA Free : Tritan body protects from kids to adults!
- Easy to Use : Portable & Chargeable & Unbreakable

Amazing ③ in ① Effect

- ① Drinking H₂ Water : Human Body Anti-aging Effect
- ② With Mineral Cartridge(Optional) : Hydrogen Mineral Effect !
- ③ With H₂ Water Mist : Facial Anti-aging Effect!

COMPOSITION

- Main Unit, USB Cable, Applicator, User Manual, Package Box, Mineral Cartridge(optional), Hydrogen Water Mist(Optional)

Nutricare



OEM•ODM

KEY TO SUCCESS



SUCCESS



MARKET RESEARCH

PLANNING

MANAGEMENT



Nutricare

Company Info



"Always,
Your Best Solution"



Nutricare Co., Ltd.

11F, 88, Yulgok-ro, Jongno-gu, Seoul 03131 Korea

TEL. +82.2.3672.1841

FAX. +82.2.3672.1844

www.nutricare.co.kr