


Your vision,
Our mission.





What is your style?

01

COVER UP
NATURAL LENS

We have more than 3,000 different and trendy lenses which are suitable for each country. These lenses can make your iris more naturally and shape of your eyes more clearly.



PEARL LENS

Glittering gold dust, silver dust and unique lens.




FANCY LENS


Various lenses of artistic design that can be used for special makeup to Halloween party.



IRIS LENS


Beauty lens for people with iris color problems.






DCP Method (Double Coating Printing)

ICK's DCP method prints a pigment between contact lens material, and keeps the pigment from coming in direct contact with the eye, providing comfortable and safe color contact lenses.



Oxygen Permeability

ICK is a contact lens manufacturer that continuously strives to produce comfortable and healthy-looking contact lenses. The silicone hydrogel contact lenses developed by the company R&D center offers high oxygen permeability, excellent protein non-adherence, consistent comfort and moisture for your eyes all day long without irritating deposits.



True-eye Edge Design

In order to provide a comfortable fit, the part of ICK contact lenses that comes in contact with the cornea has an ergonomic true eye edge design that fits comfortably. The circle line between the lens material enables breathability causing no side effects due to lack of oxygen.

Optic Zone : 8.0mm ~ 9.5mm

Printing Zone : 13.3mm ~ 14.0mm

Base Curve : 8.6mm ~ 8.8mm

Diameter : 14.0mm ~ 15.0mm

Center Thickness : 0.05~0.12mm

PACKAGE STYLE



OPTIMUS
52*26*105(mm)



BASIC
45*88*15(mm)



CUSTOM



MINI BLISTER PACK



SEAL LINER



LABEL

Lens Type			D	W	M	Y
CLEAR	Myopia	Polymacon	✓	✓	✓	✓
		Silicon Hydrogel	✓	✓	✓	
	Hyperopia	Polymacon			✓	✓
		Silicon Hydrogel				
	Astigmatism	Polymacon			✓	✓
		Silicon Hydrogel			✓	
COLOR	Myopia	Polymacon	✓	✓	✓	✓
		Silicon Hydrogel	✓	✓	✓	
	Hyperopia	Polymacon			✓	✓
		Silicon Hydrogel				
	Astigmatism	Polymacon			✓	✓
		Silicon Hydrogel				
PEARL	Myopia	Polymacon	✓	✓	✓	
		Silicon Hydrogel				
FANCY	Spherical Type				✓	✓
	Prism Type				✓	✓
IRIS	Polymacon					✓

Myopia : 0.00~ -5.00 D : -0.25 Step, -5.00 ~ -12.00 D : -0.50 Step (RX : -8.50 ~ -12.00) / DIA 14.0~15.0 / BCR : 8.6mm

Hyperopia : 0.00~ +5.00 D : +0.25 Step, +5.00 ~ +6.00 D : +0.50 Step / DIA 14.0~15.0 / BCR : 8.6mm

Toric : ① Sphere Power : 0.00~ -5.00 D : -0.25 Step, -5.00 ~ -12.00 D : -0.50 Step
0.00~ +5.00 D : +0.25 Step, +5.00 ~ +6.00 D : +0.50 Step


② Cylinder Power : -0.50~-4.50 D : -0.25 Step / Axis : 10~180 (10 degree) / Diameter : 14.0mm / BCR:8.6mm

03 ICK's unique manufacturing process!

All these processes take place within 45 days.

within 20 to 25 days.

within 20 to 25 days.





ICK contact lenses create a beautiful world by striving to achieve zero customer complaints through a systematic quality management and 100% total inspection.



2F, Hyundai Town, 99, Gwangdeokdong-ro, Danwon-gu, Ansan-si,
Gyeonggi-do, Republic of Korea (15462)



36, Madogongdan-ro 2-gil, Mado-Myeon, Hwaseong-si,
Gyeonggi-do, 18542, Rep. of KOREA



+82-31-408-8901



sales@icklens.com



www.icklens.com