



## 

WATER SOFTENER

















ung WATER SOFTENER

UNG WATER SOFTENER

Size(Diameter6.3xLength15.3cm) Specifications

Weight(350g)

Color(Silver, Gold, Black, Pearl)

IP Rating

Manufacturer / Origin



urbo self powered generator



UNG (주)유엔지 테크놀로지 U&G Technology Co.,Ltd.

Office hours | Monday ~ Friday (9:30 ~ 17:30) Lunch hours | 12:00 ~ 13:00


(closed on weekends and holiday

www.ungth.com



## **Electronic Charming Water!**

This is the smallest **UNG** Water Softener using electrostatic field. Improve the quality of incoming water and provides a high quality water environment,






## Scientifically designed home skin manager **UNG** Water Softener!

Change of Water, Skin · Hair PLUS CARE

FEATURES





15% softer water than conventional tap water when using this product.



effect

Prevents dermatitis and skin diseases by reducing 66.6% of bacteria.



cleaning

Cleans up dirt contaminated in pores, such as fine dust and vellow dust!



moisturizing

Water molecules with weakened hydrogen bonds penetrate deep into pores.



Skin keratin control

Improves skin keratin by skin moisturizing



Scalp control

Effective removal of residual detergents improves scalp and hair



- ☐ The Skin is dry, rough, and itchy. ☐ The hair is stiff and lacking in luster.
- ☐ The soap is not removed cleanly on the skin.

We recommend

**UNG WATER SOFTENER to** 

- ☐ Difficult to remove stains on tiles, glass, and bathtubs.
- The faucet is covered with white lime.
- ☐ After brushing teeth, toothpaste seems to have left in mouth, feeling dry.
- ☐ Even if you wash your feet, smell from the feet.
- ☐ To bathe pets more cleanly.

## How to install **UNG** Water Softener

powered

Economical

Convenient.



Remove the screw of the shower hose connected to the shower faucet by a spanner.



Connect the inlet of the product to the separated shower faucet.



Connect the shower hose to the oulet of the product.



Make sure the water is not leaking and connected properly.

Korea Conformity Laboratories